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Elektrostymulator TENS EMS AS-2016

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User manual

Electrostimulator TENS EMS AS-2016

The AS-2016 TENS EMS electrostimulator is a **portable therapeutic, rehabilitation, and relaxation device** that **naturally relieves pain, inflammation, and swelling.**

Please read this user manual before use.

Dear Sir or Madam, thank you very much for purchasing our product.

This manual has been developed for the AS-2016 electrostimulator. It describes the electrostimulator, its application, menu, functions, technical specifications, how to use it, troubleshooting, etc. The manual also describes safety procedures designed to protect both the user and the device itself. Detailed information can be found in the relevant chapters.

Please read this manual carefully before using the device.

Failure to follow the instructions may result in damage to the equipment or personal injury. The manufacturer, importer, or distributor are not responsible for the safety, reliability, performance, personal injury, or damage to the device resulting from the user's failure to follow the instructions in this manual. Warranty service does not cover such defects. This product is designed for repeated use.

Electrostimulation is most commonly used in orthopedics, neurology, and rheumatology, as well as in medical rehabilitation for various degenerative diseases of the joints, bones, and muscles. Muscle electrostimulation is an electrotherapeutic treatment, also known as electrogymnastics, which aims to elicit a motor response in a muscle or group of muscles.

The purpose of electrostimulation is to regenerate muscles or muscle groups, increase their strength and endurance, increase or decrease muscle and joint tension, improve circulation, protect against atrophy by improving local circulation and increasing temperature, reduce swelling, and reduce spasticity. Electrostimulation involves pulses applied using point electrodes. The aim of the treatment is to stimulate the muscles to contract (the muscles contract and relax in rhythm with the intermittent current).

Electrostimulator user manual

1. The device may need to be charged for about an hour when the battery level is low.
2. Connect the cable to the output port at the bottom of the device (A1, A2, B1, B2).
Each port can be used individually, allowing you to place 2, 4, 6, or 8 electrodes on your body.
3. Connect the electrodes to the cable by snapping them into place.
4. Remove the protective film and place the electrodes on your body, making sure they are on your skin and do not overlap.
5. Turn on the device using the slider (ON/OFF).
6. When turned on, the device will start in program 1.
7. Select one of the 24 programs on the display. After selecting a program, gradually increase the intensity by pressing “+” or decrease it by pressing “-”. The power adjustment buttons are located at the bottom of the device. The power of ports A and B is adjusted separately. The A/B button switches the control between ports.
8. To set the time, press the clock icon on the display and select a time between 10 and 60 minutes. The default setting is 20 minutes. **IMPORTANT: The device does not work without electrodes connected to the body.**
9. After the treatment, turn off the device with the slider (OFF) and then disconnect the electrodes.
10. Protect the electrodes with the original protective film or place them on a double-sided pad.
11. The “key” button locks all functions and program selection.
12. The “pause” button pauses the device while retaining the previous settings.

Description of device functions

1. Time setting: Tap the clock icon to set the time: 10, 20, 30, 40, 50, or 60 minutes.
2. Output A or B change field: Tap the B icon to change the channel from A to B, or tap the A icon to change the channel from B to A.
3. Screen lock icon: Press the key icon once to lock the touch screen; press again to unlock it. When the screen is locked, the lock icon flashes.
4. Output A
5. USB port – for charging the battery
6. Output B
7. On/off switch – slide the switch left or right to turn the device on/off
8. Battery level indicator
9. “+” button – increases the intensity (1–20 levels)
10. “-” button – decreases intensity (1–20 levels)
11. Program selection – select one of 24 therapy programs on the touch screen
12. Pause – pauses all device operations
13. Output A/B change button

Fig. description of device functions

Mode	Name	Description	Function
1.	Tuina	Tuina simulation is a practical body treatment that opens up the flow of Qi through the meridians.	Particularly effective for joint pain, sciatica, muscle spasms, and pain in the back, neck, and shoulders. It also helps with insomnia, constipation, headaches, and tension.
2.	Acupuncture	Simulated acupuncture feels like the real thing, without the use of needles, using pulses to various depths of puncture.	By stimulating specific acupuncture points, it corrects disturbances in the flow of Qi through the meridians.
3.	Tapping	It simulates impacts that massage the muscles, giving a feeling of invigoration and revitalization.	Particularly effective in relaxing muscles.
4.	Guasha	Low-frequency currents. They provide a scratching and kneading sensation. They release unhealthy elements from damaged areas, promoting circulation and healing.	Effective in acute and chronic internal diseases. Provides immediate relief from pain, stiffness, coughing, nausea, etc.
5.	Cupping	Bubble simulation, constant stimulation to tense the muscles.	Relieves back pain, neck pain, muscle stiffness, migraines, fatigue, rheumatism, and cellulite.
6.	Shiatsu	A "finger pressure" massage that restores the body's energy balance, provides deep relaxation, and reduces muscle tension.	It helps with headaches, back pain, digestive problems, sinus issues, and skin problems.
7.	Pressure	Manipulation of muscles and connective tissue using various techniques.	It improves functions, promotes healing, reduces muscle tension, relaxes, and improves well-being.
8.	Back Pain	Low-frequency vibration promotes circulation and relieves pain.	Especially recommended for the lower back.
9.	Deep Tissue Massage	Deep tissue massage/Swedish massage. Strong pressure relieves chronic muscle tension.	Provides relief from pain and stiffness associated with sports injuries, strain, degeneration, muscle spasms, and carpal tunnel syndrome.
10.	Body Building	Electrical impulses cause muscle contractions, which leads to their strengthening.	It helps build muscle, relax spasms, increase range of motion, and aid recovery after injuries and surgeries.
11.	Weight Loss	A combination of acupuncture and impulses.	It helps you lose weight and strengthen and tone the muscles of your abdomen, back, thighs, and arms.
12.	Combination	Automatic combination of modes for relaxation.	Complete relaxation of the body after training, exercise, or a hard day's work.

13. Shoulder 	14. Neck 	15. Back 	16. Elbow 	17. Hip 	18. Ankle 
19. Foot 	20. Hand 	21. Knee 	22. Wrist 	23. Ear 	24. Relax 

Operational safety

- ⊗ Never attach the electrodes to your skin while the device is turned on – there is a risk of electric shock.
- ⊗ Dispose of the device, accessories, batteries, and packaging in accordance with local regulations.
- ⊗ If the device malfunctions or causes discomfort, stop using it immediately.
- ⊗ Do not use the device in a damp environment, while bathing, driving, or sleeping.
- ⊗ Service and maintenance may only be performed by qualified personnel. The user must not perform maintenance themselves.
- ⊗ Do not use the electrostimulator simultaneously with devices not described in this manual.
- ⊗ Do not connect the electrodes to metal devices.
- ⊗ Keep the device away from corrosive substances, explosives, high temperatures, and moisture.
- ⊗ If the electrostimulator or accessories get wet, stop using them immediately.
- ⊗ After moving the device from a cold and humid environment to a warm one, wait before turning it on.
- ⊗ Do not disinfect the device with hot steam or high pressure.
- ⊗ Do not immerse the electrostimulator in liquids – clean it with a soft cloth soaked in medical alcohol.
- ⊗ Do not use the device in a flammable environment risk of explosion.
- ⊗ Store the device in a dry place – excessive moisture may damage it.
- ⊗ Do not press the main panel with sharp objects.

Warnings and contraindications

- ⊗ Heart disease without prior consultation with a physician.
- ⊗ Persons with pacemakers or other implanted defibrillation devices.
- ⊗ Diseases such as epilepsy, arrhythmia.
- ⊗ High fever, chronic infectious diseases.
- ⊗ Early period after muscle injury (7–10 days).
- ⊗ Abnormal blood pressure.
- ⊗ Menstruation, pregnancy, or postpartum period.
- ⊗ Breastfeeding.
- ⊗ Abdominal or inguinal hernia.
- ⊗ Benign and malignant tumors.
- ⊗ On cuts / skin diseases / skin lesions / edema / phlebitis / thrombophlebitis / varicose veins.
- ⊗ If the body is wet after bathing, showering, or as a result of sweating.
- ⊗ Tuberculosis.
- ⊗ Presence of metal foreign bodies in tissues.
- ⊗ Blood vessel insufficiency.
- ⊗ Inflammation of the skin and soft tissues.
- ⊗ Eczema (dermatological lesions).
- ⊗ Electrodes must not be placed in the area of the carotid sinus, neck, larynx, mouth, and throat, or on the lips, as this may cause severe spasms of the larynx and throat muscles; these spasms may be severe enough to close the airways or cause breathing difficulties.
- ⊗ Do not place electrodes above or on the eyes.
- ⊗ Do not stimulate the chest in the area of the heart, as the introduction of electrical current into the heart may cause cardiac arrhythmia.
- ⊗ Stimulation must not be applied through the brain.

Additional comments


- ⊗ The long-term effects of chronic electrical stimulation are unknown.
- ⊗ Do not use the device on children, infants or newborns.
- ⊗ The product is intended for adults.
- ⊗ If you are unsure whether you can use the electrostimulator, consult your doctor.

Indications for use

- ⊗ body shaping
- ⊗ weight loss (fat reduction)
- ⊗ increasing muscle endurance
- ⊗ relaxation and unwinding
- ⊗ warm-up before training
- ⊗ muscle development (abdominal muscles, deltoids, biceps, triceps, forearm muscles, glutes, quadriceps, calf muscles)
- ⊗ prevention or inhibition of muscle spasms
- ⊗ re-education of lost muscle functions
- ⊗ maintaining or increasing the range of motion in the joints
- ⊗ rehabilitation after injuries
- ⊗ urinary incontinence-stimulation of the peroneal nerve
- ⊗ back pain (cervical, shoulder, thoracic, lumbar)
- ⊗ Torticollis (stiffness/spasm of the neck)
- ⊗ Improvement of local blood circulation
- ⊗ Pain after mechanical injuries
- ⊗ Degenerative conditions
- ⊗ Neuralgia
- ⊗ Rheumatism
- ⊗ Migraines
- ⊗ Tennis elbow syndrome
- ⊗ Elbow pain
- ⊗ Knee pain
- ⊗ Foot drop

2 Cables


1 Cables

TENS


Chronic cervical spine pain

⊖ PROGRAM: 1,2,3,4,5,6,7,9,13

Wrist pain

⊖ PROGRAM: 1,2,3,4,16,22


Ankle pain

⊖ PROGRAM: 1,2,3,5,18

Relieving stomach ache

⊖ PROGRAM: 1,7


Elbow joint pain

⊖ PROGRAM: 1,2,3,4,7,9,16,22


Knee pain

⊖ Osteoarthritis of the knee
⊖ PROGRAM: 1,2,3,4,7,9,10,21


Shoulder pain

⊖ Compression of the spinal cord in the cervical region
⊖ PROGRAM: 1,2,3,4,5,6,7,13,22


Excessive tension in the paraspinal muscles

⊖ Thoracic section
⊖ PROGRAM: 1,2,3,4,7,9,10,21


Sciatica

⊖ PROGRAM: 1,2,3,4,5,8,17


Shoulder pain

⊖ PROGRAM: 2,3,4,13


Pain or tension in the calf

PROGRAM: 1,2,3,4,5,6,7,18


Ankle pain

PROGRAM: 1,2,3,4,18

Shoulder blade pain


PROGRAM: 1,2,3,4,5,6,9,15

Lumbosacral pain

Tension in the lumbar spine muscles without pain radiating to the lower limbs


PROGRAM: 1,2,3,4,5,6,7,8,9,15

Tension in the paraspinal muscles


Cervical section

PROGRAM: 1,2,3,4,5,6,7,13,15


Hip pain

PROGRAM: 1,2,3,4,5,6,7,9,17


Urinary incontinence

- ⊖ Electrostimulation of the tibial nerve (short and long) (peroneal)
- ⊖ PROGRAM: 10


Gluteal muscles

- ⊖ PROGRAM: 10,11


Abdominal muscles

- ⊖ PROGRAM: 10,11


Falling foot

- ⊖ Damage to the nerve or atrophy of the tibialis anterior muscle
- ⊖ PROGRAM: 10


Quadriceps muscle of the thigh

- ⊖ PROGRAM: 10,11

Deltoid muscle

- ⊖ PROGRAM: 10,11


Arm muscles

PROGRAM: 10,11


Forearm muscles

PROGRAM: 10,11

Calf muscles

PROGRAM: 10,11

Biceps muscle

PROGRAM: 10,11

ATTENTION!

Electrotherapy does not replace medical treatment.

Before using the electrostimulator, it is recommended to consult a doctor or physiotherapist, especially in the case of chronic diseases, pain of unknown origin, after injuries or surgical procedures.

The user uses the device at their own risk and is obliged to follow the safety rules and recommendations contained in the user manual. The manufacturer, importer or distributor are not liable for the consequences of improper use of the device or use that is not in accordance with the instructions.

How to care for electrodes to extend their service life

- ⊗ Make sure the device is turned off.
- ⊗ If necessary, remove hair from the application site.
- ⊗ Degrease and dry the skin (this extends the life of the electrodes).
- ⊗ Connect the electrodes to the device and remove the protective film (keep it for reuse; if damaged, use the double-sided pads from the kit).
- ⊗ Stick the electrodes to the skin and press lightly.
- ⊗ Read the device manual before use.
- ⊗ After the treatment, slowly peel off the electrodes, starting from the corner, and stick them back onto the protective film.
- ⊗ Store the electrodes in their packaging in a dry place.
- ⊗ Replace the electrodes when they stick poorly or the gel is dry or dirty.
- ⊗ Do not stick exposed electrodes together as this may damage them.
- ⊗ The electrodes are not intended for disinfection; each user should have their own.

Technical specifications of the device

- 2 stimulation techniques: EMS and TENS
- 24 programs
- 20 intensity levels / independent power adjustment for zones A and B
- Up to 8 electrodes can be connected at the same time
- Large, touch-sensitive, easy-to-read, backlit LCD display
- Operating time adjustment (10–60 minutes)
- Treatment time preview
- Touch screen lock function
- Battery power supply: 300 mAh lithium-ion battery
- Frequency 1–120 Hz, device power from 6 to 140 mA
- Charging via USB-C port (240V power supply)
- Electrode connection standard: SNAP 3.5 mm
- Certificates: CE marked device

Problem solving

Weak stimulation

- ⊗ Make sure that the electrodes do not overlap.
- ⊗ Check that the electrodes and cable are properly connected.
- ⊗ Check the correct placement of the electrodes.
- ⊗ Make sure the electrodes adhere well to the skin and are not dry; replace them with new ones if necessary.
- ⊗ Check that the stimulation power level is not set too low.

No stimulation

- ⊗ Check that the electrodes and cable are connected correctly.
- ⊗ Make sure that the electrodes are firmly attached to the skin.

One electrode is stronger than the other.

- ⊗ This is normal; different parts of the body react differently.
- ⊗ Make sure the electrodes are not dry and adhere well to the skin.

The display and therapy mode are working, but there is no stimulation

- ⊗ Check that the electrodes are positioned correctly and adhere well to the skin (the device may switch off after approx. 10 seconds if an electrode falls off or if not all of the two or four electrodes are attached to the body).
- ⊗ Check the power level and adjust the intensity so that you feel a gentle stimulation.
- ⊗ Make sure that the electrodes do not overlap.

When switched on, neither the display nor the therapy modes light up

- ⊗ Check that the device is charged.

The skin becomes reddened.

- ⊗ Check the treatment time; do not exceed 10–15 minutes.
- ⊗ Make sure that the electrodes are not dry or dirty; replace them with new ones if necessary.
- ⊗ Check that the electrodes adhere well to the skin.

Gel electrodes do not stick to the skin even after cleaning and moistening.

- ⊗ Replace the adhesive gel pads.

Accessories

- Ø 240V power supply
- Ø Electrode connecting cable – 4 pcs.
- Ø Power supply connection cable via USB-C
- Ø Double-sided pad
- Ø Stylus
- Ø Cardboard packaging
- Ø Hydrogel electrodes for the body
- Ø Electrodes – gloves – optional
- Ø Electrodes – socks – optional
- Ø Knee electrodes – optional
- Ø Electrodes – flip-flops – optional
- Ø User manual

The **AS-2016** electrostimulator is a product that meets all safety standards. It has **two CE certificates (EMC and LVD)**. The **EMC** certificate complies with Directive 2014/30/EU of the European Parliament and of the Council of 26 February 2014 on electromagnetic compatibility. The **LVD** certificate complies with Directive **2014/35/EU** of the European Parliament and of the Council of 26 February 2014 – Low Voltage Directive.

This symbol indicates that this product must not be disposed of with other household waste. The user is responsible for delivering the used device to a collection point for the disposal of electrical and electronic equipment. Separate collection and recycling of devices helps to protect the environment and human health. For more information on the disposal of devices, please contact your local authorities, disposal companies, or the place where you purchased this product.

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